

PRODUCT SPECIFICATIONS

PART NUMBER:	A123
LOT NUMBER:	D2081
DESCRIPTION:	Anti-Human IgG-Fc Fragment IgG Fraction
HOST ANIMAL:	Goat
LOT NUMBER:	D2081
TOTAL PROTEIN:	55mg/ml
SPECIFIC ANTIBODY TITER:	11mg/ml
BUFFER:	Phosphate Buffered Saline
PRESERVATIVE:	0.1% Sodium Azide (NaN ₃)
ANTIGEN USED:	Normal IgG, Fc Fragment (papain digest)

This product is being sold for research purposes only. It is not to be used in humans or animals. There is no express or implied warranty. No liability is assumed for fitness of purpose and merchantability, or direct and consequential damages. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations. 12/04