

ANTIBODY SPECIFICATION SHEET

PRODUCT:	Monoclonal THC
CATALOG NUMBER:	AG309
LOT NUMBER:	J2302
FORMAT:	Protein A Purified Material
HOST:	Murine
ANTIGEN:	8-THC-BSA
ANTIBODY SPECIFICITY:	THC and Metabolite
CONCENTRATION:	4.16 mg/mL by A280 nm
ISOTYPE:	IgG-1
BUFFER:	0.05 M Phosphate 0.85% Sodium Chloride, pH 7.15
PRESERVATIVE:	0.05% Sodium Azide
STORAGE:	Short Term -- Refrigerated Long Term -- Frozen
CROSS REACTIVITY:	No cross reactivity with: <ul style="list-style-type: none">• Opiates• Cocaine Metabolite• Amphetamines• Phencyclidine

For research and manufacturing use only ORDERING INFORMATION

Orders may be placed by phone, fax, or mail.
Please provide Part Number, Description, Quantity Desired, and Purchase Order Number.
THANK YOU. YOUR BUSINESS IS APPRECIATED.