

**ANTIBODY SPECIFICATION SHEET**

**PRODUCT:** Monoclonal Beta HCG--CTP  
**CATALOG NUMBER:** B01-96-01M-P  
**LOT NUMBER:** K10  
**FORMAT:** Purified Material  
**HOST:** Murine  
**ANTIGEN:** Beta HCG-CTP Thiogobulin  
**ANTIBODY SPECIFICITY:** Beta HCG-CTP (Carboxy-Terminal-Peptide)  
**CONCENTRATION:** 3.91 mg/ml  
by A280nm  
**Ka:** Half Displacement with:  
Whole HCG @ 15 ng/ml  
**ISOTYPE:** IgG-2a  
**BUFFER:** 0.05M Phosphate  
0.85% Sodium Chloride, pH 7.15  
**PRESERVATIVE:** 0.05% Sodium Azide  
**STORAGE:** Short Term -- Refrigerated  
Long Term -- Frozen  
**CROSS REACTIVITY:** LH =0.1%

***MANUFACTURED BY BIOSTRIDE, INC.  
FOR RESEARCH AND MANUFACTURING USE ONLY***

**ORDERING INFORMATION**

Orders may be placed by phone, fax, or mail.  
Please provide Part Number, Description, Quantity Desired, and Purchase Order Number.  
THANK YOU. YOUR BUSINESS IS APPRECIATED.