

ANTIBODY SPECIFICATION SHEET

PRODUCT: Monoclonal Beta HCG
CATALOG NUMBER: B01-96-08M-P
LOT NUMBER: E21
FORMAT: Protein A Purified Material
HOST: Murine
ANTIGEN: HCG
ANTIBODY SPECIFICITY: Beta HCG
CONCENTRATION: 5.93 mg/ml
by A280nm
KA: 1X 10¹¹ Liters/Mole
ISOTYPE: IgG-1
BUFFER: 0.05M Phosphate
0.85% Sodium Chloride, pH 7.15
PRESERVATIVE: 0.05% Sodium Azide
STORAGE: Short Term -- Refrigerated
Long Term -- Frozen
CROSS REACTIVITY: Beta HCG = 100%
HCG = 100%
Alpha HCG = 1.1%
LH = 0.4%
FSH = 0.3%
TSH = 0.08%

MANUFACTURED BY BIOSTRIDE, INC.
For research and manufacturing use only

ORDERING INFORMATION

Orders may be placed by phone, fax, or mail.
Please provide Part Number, Description, Quantity Desired, and Purchase Order Number.
THANK YOU. YOUR BUSINESS IS APPRECIATED.