

ANTIBODY SPECIFICATION SHEET

PRODUCT: Monoclonal Beta LH
CATALOG NUMBER: B01-96-42M-P
LOT NUMBER: R22
FORMAT: Protein A Purified
HOST: Murine
ANTIGEN: LH
ANTIBODY SPECIFICITY: Human Urinary LH
CONCENTRATION: 2.65 mg/ml
by A280nm
KA: 1.8 X 10¹¹ Liters/Mole LH
ISOTYPE: IgG-1
BUFFER: 0.05M Phosphate
0.85% Sodium Chloride, pH 7.15
PRESERVATIVE: 0.05% Sodium Azide
STORAGE: Short Term -- Refrigerated
Long Term -- Frozen

CROSS REACTIVITY:

Beta LH	= 55%
FSH	=0.9%
TSH	=18.4%
Alpha HCG	=4%
HCG	=57%

ORDERING INFORMATION

Orders may be placed by phone, fax, or mail.
Please provide Part Number, Description, Quantity Desired, and Purchase Order Number.
THANK YOU. YOUR BUSINESS IS APPRECIATED.