

ANTIBODY SPECIFICATION SHEET

PRODUCT:	Monoclonal Morphine
CATALOG NUMBER:	B01-99-02M-P
LOT NUMBER:	P04
FORMAT:	Protein A Purified Material
HOST:	Murine
ANTIGEN:	Morphine-BSA
ANTIBODY SPECIFICITY:	Morphine
CONCENTRATION:	5.47 mg/ml by A280nm
KA:	Half Displacement with Morphine @ 0.4ng/ml
ISOTYPE:	IgG-1
BUFFER:	0.05M Phosphate 0.85% Sodium Chloride, pH 7.15
PRESERVATIVE:	0.05% Sodium Azide
STORAGE:	Short Term -- Refrigerated Long Term -- Frozen

***MANUFACTURED BY BIOSTRIDE, INC.
For research and manufacturing use only***

ORDERING INFORMATION

Orders may be placed by phone, fax, or mail.

Please provide Part Number, Description, Quantity Desired, and Purchase Order Number.
THANK YOU. YOUR BUSINESS IS APPRECIATED.



EastCoast Bio, Inc.
211 Wells Street
P.O. Box 489
North Berwick, ME 03906

P 207.676.7639
F 207.676.7658
E info@eastcoastbio.com
www.eastcoastbio.com