

**ANTIBODY SPECIFICATION SHEET**

|                        |  |
|------------------------|--|
| <b>PRODUCT:</b>        | Monoclonal Hepatitis B Surface Antigen                   |
| <b>CATALOG NUMBER:</b> | B04-95-33M-P   |
| <b>LOT NUMBER:</b>     | F12  |
| <b>FORMAT:</b>         | Purified Material  |
| <b>HOST:</b>           | Murine   |
| <b>ANTIGEN:</b>        | HbsAg  |
| <b>CONCENTRATION:</b>  | 5.60 mg/mL<br>by A280 nm (extenction coefficient = 1.08) |
| <b>Ka:</b>             | 1 x 10 <sup>10</sup> L/M                                 |
| <b>ISOTYPE:</b>        | IgM  |
| <b>BUFFER:</b>         | 0.05 M Phosphate<br>0.85% Sodium Chloride, pH 7.15       |
| <b>PRESERVATIVE:</b>   | 0.05% Sodium Azide                                       |
| <b>STORAGE:</b>        | Short Term -- Refrigerated<br>Long Term -- Frozen        |

***MANUFACTURED BY BIOSTRIDE, INC.***  
**for research and manufacturing use only**  
**ORDERING INFORMATION**

Orders may be placed by phone, fax, or mail.  
Please provide Part Number, Description, Quantity Desired, and Purchase Order Number.  
THANK YOU. YOUR BUSINESS IS APPRECIATED.