

ANTIBODY SPECIFICATION SHEET

PRODUCT: Monoclonal Methamphetamine
CATALOG NUMBER: C03-99-31M-P
LOT NUMBER: R11
FORMAT: Protein A Purified Material
HOST: Murine
ANTIGEN: NA
ANTIBODY SPECIFICITY: Methamphetamine
PURITY: >98%
CONCENTRATION: 5.2 mg/ml
by A280nm
ISOTYPE: IgG2a,k
BUFFER: PBS, pH 7.2
PRESERVATIVE: 0.05% Sodium Azide
STORAGE: Short Term -- Refrigerated
Long Term -- Frozen

ORDERING INFORMATION**For research and manufacturing use only**

Orders may be placed by phone, fax, or mail.
Please provide Part Number, Description, Quantity Desired, and Purchase Order Number.
THANK YOU. YOUR BUSINESS IS APPRECIATED.