

ANTIBODY SPECIFICATION SHEET

PRODUCT: Monoclonal Amphetamine
CATALOG NUMBER: C03-99-72M-P
LOT NUMBER: J24
FORMAT: Protein A Purified Material
HOST: Murine
ANTIGEN: N/A
ANTIBODY SPECIFICITY: Amphetamine
CONCENTRATION: 4.9 mg/ml
by A280nm
ISOTYPE: IgG2b, k
BUFFER: Phosphate Buffered Saline, pH 7.2
PRESERVATIVE: 0.05% Sodium Azide
STORAGE: Short Term -- Refrigerated
Long Term -- Frozen

MANUFACTURED BY BIOSTRIDE, INC.
for research and manufacturing use only
ORDERING INFORMATION

Orders may be placed by phone, fax, or mail.
Please provide Part Number, Description, Quantity Desired, and Purchase Order Number.
THANK YOU. YOUR BUSINESS IS APPRECIATED.