

POLYCLONAL ANTIBODY DATA SHEET

PRODUCT:	Rabbit polyclonal anti <i>E.coli</i> O157
CATALOG NUMBER:	C101C-G2611
CONCENTRATION:	1.07 mg/ml
PRESENTATION:	Affinity purified antibody in PBS with 0.05% sodium azide
STORAGE:	Short: Refrigerate Long: Freeze Avoid repeated freezing and thawing.

This product is being sold for research or further manufacturing use only. It is not to be used in humans or animals. There is no express or implied warranty. No liability is assumed for fitness of purpose and merchantability, or direct and consequential damages. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.8/02