

CMV ICP36 Monoclonal Antibody

Catalog Number	CA006
Lot Number	G0729146
Expiration Date	July 2010
Target Protein	Cytomegalovirus ICP36 DNA Binding Protein
Target Protein Also Known As	ppUL44, MDBP, pp52, DB52
Clone Designation	10D8
Species/Isotype	Mouse IgG1(k)
Purity	Affinity Purified, >95% IgG by SDS-PAGE
Supplied In	Phosphate Buffered Saline, pH 7.4 (no sodium azide)
Protein Concentration	1.0 mg/mL
Immunofluorescence Titer	1+ @ 1:12,800
ELISA Titer	0.100 @ >1:102,400
Western Blot Titer	>1:80,000
Reference and/or Additional Reading	J. Virology (1995) 69(4):2047

This product is intended for research and manufacturing uses only. It is not a diagnostic device. It is supplied frozen on dry ice. Upon receipt, store at -20°C. Product degradation will result from multiple freeze/thaw cycles. It is suggested that the antibody be stored in use size aliquots and thawed just prior to use. Storage at 2-8°C for up to 2 weeks is acceptable for unused portions.

Import/Export Notice: This product contains no serum or plasma of animal origin. This product conforms to the specifications set forth in PS-106 and is suitable for release.