

ANTIBODY SPECIFICATION SHEET

PRODUCT:	Monoclonal anti-Propoxyphene
CATALOG NUMBER:	CB302
LOT NUMBER:	P2161
HOST:	Murine
ANTIBODY SPECIFICITY:	Propoxyphene
ISOTYPE:	IgG
FORMAT:	Protein A Purified Material
CONCENTRATION:	3.45 mg/mL
BUFFER:	0.05M KPhos, 0.85% NaCl, pH 7.15
PRESERVATIVE:	0.05% Sodium Azide
STORAGE:	Short Term Refrigerated, 2-8°C Long Term Frozen, -20°C

For Research or Further Manufacturing Use Only
Not for Use in Diagnostic Procedures

ORDERING INFORMATION

Orders may be placed by phone, fax, or mail.
Please provide Part Number, Description, Quantity Desired, and Purchase Order Number.
THANK YOU. YOUR BUSINESS IS APPRECIATED.



EastCoast Bio, Inc.
211 Wells Street
P.O. Box 489
North Berwick, ME 03906

P 207.676.7639
F 207.676.7658
E info@eastcoastbio.com
www.eastcoastbio.com