

**POLYCLONAL ANTIBODY DATA SHEET**

<b>PRODUCT:</b>	Anti-Mouse IgG (H & L)
<b>CATALOG NUMBER:</b>	F101
<b>LOT NUMBER:</b>	Q2212
<b>HOST:</b>	Goat
<b>FORMAT:</b>	Affinity Purified
<b>CONCENTRATION:</b>	5.0 mg/ml
<b>PURITY:</b>	>95% by SDS-PAGE
<b>BUFFER:</b>	PBS, pH 7.2
<b>PRESERVATIVE:</b>	0.1% sodium azide
<b>STORAGE:</b>	Short term: 4°C Long term: -20°C
<b>SPECIFICATIONS:</b>	This product reacts with all subclasses of IgG and with all classes of mouse immunoglobulins by immunodiffusion and IEP techniques due to the presence of anti-light chains. This product is suitable for ELISA and Blotting applications.
<b>EXPIRATION DATE:</b>	August 2012

*This product is intended for research and manufacturing uses only. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.*

**ORDERING INFORMATION**

Orders may be placed by phone, fax, or mail.  
Please provide Part Number, Description, Quantity Desired, and Purchase Order Number.  
THANK YOU. YOUR BUSINESS IS APPRECIATED.



**EastCoast Bio, Inc.**  
211 Wells Street  
P.O. Box 489  
North Berwick, ME 03906

P 207.676.7639  
F 207.676.7658  
E [info@eastcoastbio.com](mailto:info@eastcoastbio.com)  
[www.eastcoastbio.com](http://www.eastcoastbio.com)