

HSV-1 gE Monoclonal Antibody

Catalog Number	H1A054
Lot Number	H0730107
Expiration Date	August 2010
Target Protein	Herpes Simplex Type 1 Glycoprotein E
Target Protein Also Known As	gE-1, pUs9, 80 kD Fc Receptor
Clone Designation	9H3
Species/Isotype	Mouse IgG2a(k)
Purity	Affinity Purified, >95% IgG by SDS-PAGE
Supplied In	Phosphate Buffered Saline, pH 7.4 (no sodium azide)
Protein Concentration	1.0 mg/mL
Immunofluorescence Titer	1+ @ 1:6,400
ELISA Titer	0.100 @ >1:102,400
Reference and/or Additional Reading	J. Virol. (1987) 61(7):2208

This product is intended for research and manufacturing uses only. It is not a diagnostic device. It is supplied frozen on dry ice. Upon receipt, store at -20°C. Product degradation will result from multiple freeze/thaw cycles. It is suggested that the antibody be stored in use size aliquots and thawed just prior to use. Storage at 2-8°C for up to 2 weeks is acceptable for unused portions.

Import/Export Notice: This product contains no serum or plasma of animal origin. This product conforms to the specifications set forth in PS-154 and is suitable for release.