

HSV-2 gG Monoclonal Antibody

| | |
|--|--|
| Catalog Number | H2A023 |
| Lot Number | J0833012 |
| Expiration Date | October 2011 |
| Target Protein | Herpes Simplex Type 2 Glycoprotein G |
| Target Protein Also Known As | pUS4, gG-2 |
| Clone Designation | 3D4 |
| Species/Isotype | Mouse IgG1(k) |
| Purity | Affinity Purified, >95% IgG by SDS-PAGE |
| Supplied In | Phosphate Buffered Saline, pH 7.4 (no sodium azide) |
| Protein Concentration | 1.0 mg/mL |
| Immunofluorescence Titer | 1+ @ 1:6,400 |
| ELISA Titer | 0.100 @ >1:102,400 |
| Western Blot Titer | >1:16,000 |
| Reference and/or Additional Reading | J. Virol. (1984) 50(2):547 |

This product is intended for research and manufacturing uses only. It is not a diagnostic device. It is supplied frozen on dry ice. Upon receipt, store at -20°C. Product degradation will result from multiple freeze/thaw cycles. It is suggested that the antibody be stored in use size aliquots and thawed just prior to use. Storage at 2-8°C for up to 2 weeks is acceptable for unused portions.

Import/Export Notice: This product contains no serum or plasma of animal origin. This product conforms to the specifications set forth in PS-123 and is suitable for release.