

MONOCLONAL ANTIBODY DATA SHEET

PRODUCT: RUBELLA

CATALOG NUMBER: H301

LOT NUMBER: A1681

HOST: Mouse

CLONE: C1.19.3

PRESENTATION: ASCITES FLUID SUPPLIED IN BORATE BUFFER;
0.04 M BORIS ACID, 5MM SODIUM BORATE,
0.2 M SODIUM CHLORIDE, pH + 8.0(+/-)0.2

CONCENTRATION: 4 MG / ML

ISOTYPE: IgG2ak

COMMENTS:

- SPECIFIC FOR E1
- REACTS ON NON-REDUCED ANTIGENS.
- EXHIBITS HI LEVEL HI ACTIVITY
- NO NEUTRALIZATION OF INFECTIVITY
- SEVEN DIFFERENT SUB-TYPES

PRESERVATIVE ADDED: 0.1 % SODIUM AZIDE

This product is being sold for research purposes only. It is not to be used in humans or animals. There is no express or implied warranty. No liability is assumed for fitness of purpose and merchantability, or direct and consequential damages. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.

2/95