

HHV-8 K-bZIP Mouse Monoclonal Antibody

SPECIFICATIONS

Product Description:

Mouse monoclonal antibody to K-bZIP of Human Herpes Virus Type 8. This antibody originates from ascites fluids and is purified by protein G agarose affinity chromatography.

Catalog Number:

H8A197

Lot Number:

D0728178

Expiration Date:

April 2010

Unit Size:

500 micrograms

TECHNICAL DATA

Species/Antibody Class:

Mouse IgG

Clone Designation:

F33P1

Protein Concentration:

1.0 mg/ml (A₂₈₀)

Supplied in:

Purified IgG supplied in PBS, pH 7.4 (no azide)

Purity:

Protein G agarose affinity chromatography

Functional Activity:

Reactive with SUMOylated and non-SUMOylated species of K-bZIP of Human Herpes Virus Type 8 in Immunofluorescence (IFA) and western blot assays at 10 µg/ml.

Shipping and Storage:

This product is supplied frozen on dry ice. Upon receipt, store at -20°C. Avoid multiple freeze-thaw cycles as product degradation may result.

This product is intended for research and manufacturing uses only. It is not for use in diagnostics procedures.