

**MONOCLONAL ANTIBODY SPECIFICATION SHEET**

<b>PRODUCT</b>	Monoclonal Anti-Influenza A Virus
<b>CATALOG NUMBER</b>	J320
<b>LOT NUMBER</b>	P1066
<b>FORMAT</b>	Protein A Purified
<b>PURITY</b>	>90%
<b>HOST</b>	Mouse
<b>SPECIFICITY</b>	Nucleoprotein of many Influenza A virus strains (>50 H1 – H14 isolates)
<b>CONCENTRATION</b>	~5 mg/mL
<b>ISOTYPE</b>	IgG2a
<b>SUPPLIED IN</b>	0.01M phosphate buffered saline, pH 7.2
<b>PRESERVATIVE</b>	0.1% sodium azide
<b>APPLICATIONS</b>	IFA, IHC, ELISA (Pairs with J333)
<b>STORAGE</b>	Store at -20°C Avoid repeated freezing and thawing

*This product is intended for research and manufacturing uses only. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.*

**ORDERING INFORMATION**

Orders may be placed by phone, fax, or mail.  
Please provide Part Number, Description, Quantity Desired, and Purchase Order Number.  
THANK YOU. YOUR BUSINESS IS APPRECIATED.