

MONOCLONAL ANTIBODY SPECIFICATION SHEET

PRODUCT	Monoclonal Anti-Influenza A Virus
CATALOG NUMBER	J321
LOT NUMBER	P1067
FORMAT	Protein A Purified
PURITY	>90%
HOST	Mouse
SPECIFICITY	Nucleoprotein of all Influenza A virus strains
CONCENTRATION	~5 mg/mL
ISOTYPE	IgG2a
SUPPLIED IN	0.01M phosphate buffered saline, pH 7.2
PRESERVATIVE	0.1% sodium azide
APPLICATIONS	ELISA, IFA, IHC
STORAGE	Store at -20°C Avoid repeated freezing and thawing

This product is intended for research and manufacturing uses only. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.

ORDERING INFORMATION

Orders may be placed by phone, fax, or mail.
Please provide Part Number, Description, Quantity Desired, and Purchase Order Number.
THANK YOU. YOUR BUSINESS IS APPRECIATED.