

MONOCLONAL ANTIBODY DATA SHEET

PRODUCT:	Penicillin G
CATALOG NUMBER:	K300K
LOT NUMBER:	M0961
HOST:	Mouse
PURIFICATION:	Protein A affinity chromatography from ascites
CONCENTRATION:	1.2 mg/mL in PBS
SUBCLASS:	IgG1
IMMUNOGEN:	Penicillin-BSA conjugate
SPECIFICITY:	Reacts specifically with Penicillin G but not at the beta-lactam ring
CROSSREACTIVITY:	Does <u>not</u> react with ampicillin, cloxicillin and novobiocin
APPLICATION:	ELISA Western blot
STORAGE:	Short Term 2-8°C, Long Term -20°C

This product is being sold for research or further manufacturing use only. It is not to be used in humans or animals. There is no express or implied warranty. No liability is assumed for fitness of purpose and merchantability, or direct and consequential damages. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations. 4/07