

MONOCLONAL ANTIBODY DATA SHEET

PRODUCT:	Chlamydia LPS Mab
CATALOG NUMBER:	OO34
LOT NUMBER:	K2971
HOST:	Mouse
SUBCLASS:	IgG2b. Kappa
PRESENTATION:	Protein A purified, stored in 10mM PBS, 0.15M Sodium Chloride, pH 7.4
CONCENTRATION:	4.84 mg/mL
PRESERVATIVE:	0.1% NAN3
STORAGE:	Store at 4°C SHORT TERM -20°C LONG TERM

This product is being sold for research or further manufacturing use only. It is not to be used in humans or animals. There is no express or implied warranty. No liability is assumed for fitness of purpose and merchantability, or direct and consequential damages. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations