

ANTIBODY SPECIFICATION SHEET

PRODUCT: Monoclonal Progesterone
CATALOG NUMBER: P01-92-11M-P
LOT NUMBER: R07
FORMAT: Protein A Purified Material
HOST: Murine
ANTIGEN: Progesterone-BSA (linked at the 3 position)
ANTIBODY SPECIFICITY: Progesterone
CONCENTRATION: 1.73 mg/ml
by A280nm
SENSITIVITY: 0-50 ng/ml progesterone in EIA
ISOTYPE: IgG-1
BUFFER: PBS, pH 7.15
PRESERVATIVE: 0.05% Sodium Azide
STORAGE: Short Term -- Refrigerated
Long Term -- Frozen

CROSS REACTIVITY:

Test Compound	Progesterone
1 ug cortisol	= 5 ng
10 ng estradiol	<< 5 ng
500 ng estriol	<< 5 ng
20 ng testosterone	<< 5 ng

COMMENTS: HRP and AP conjugated matched to antibody are available.
EIA wells coated with 30 ng/well antibody onto rabbit anti-mouse

MANUFACTURED BY BIOSTRIDE, INC.
For research and manufacturer use only

ORDERING INFORMATION

Orders may be placed by phone, fax, or mail.
Please provide Part Number, Description, Quantity Desired, and Purchase Order Number.
THANK YOU. YOUR BUSINESS IS APPRECIATED.