

ANTIBODY SPECIFICATION SHEET

PRODUCT: Monoclonal Cortisol
CATALOG NUMBER: P01-92-94M-P
LOT NUMBER: Q08
FORMAT: Protein A Purified Material
HOST: Murine
ANTIGEN: Cortisol-3 CMO-BgG
ANTIBODY SPECIFICITY: Hydrocortisone
CONCENTRATION: 2.91 mg/ml
by A280nm
SENSITIVITY: D₁₀ = 0.08 ng, D₅₀ = 0.5 ng, D₉₀ = 2.8 ng
ISOTYPE: IgG
BUFFER: PBS- Azide
PRESERVATIVE: 0.05% Sodium Azide
STORAGE: Short Term -- Refrigerated
Long Term -- Frozen
CROSS REACTIVITY: Corticosterone 45%
Progesterone 9%
Deoxycortisol <2%
Dexamethazone <2%
Less than 0.01% to: Estrone (E1)
Estriol (E3)
Testosterone

***MANUFACTURED BY BIOSTRIDE, INC.
FOR RESEARCH AND MANUFACTURING USE ONLY***

ORDERING INFORMATION

Orders may be placed by phone, fax, or mail.
Please provide Part Number, Description, Quantity Desired, and Purchase Order Number.
THANK YOU. YOUR BUSINESS IS APPRECIATED.