

ANTIBODY SPECIFICATION SHEET

PRODUCT:	Monoclonal Atrazine
CATALOG NUMBER:	P01-93-41M-P
LOT NUMBER:	M39
FORMAT:	Protein A Purified
HOST:	Murine
ANTIGEN:	Atrazine- BSA
ANTIBODY SPECIFICITY:	Triazines
CONCENTRATION:	2.51 mg/ml by A280nm
ISOTYPE:	IgG-1
BUFFER:	0.05M Phosphate 0.85% Sodium Chloride, pH 7.15
PRESERVATIVE:	0.05% Sodium Azide
STORAGE:	Short Term -- Refrigerated Long Term -- Frozen

MANUFACTURED BY BIOSTRIDE, INC.

This product is intended for research and manufacturing uses only. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.

ORDERING INFORMATION

Orders may be placed by phone, fax, or mail.
Please provide Part Number, Description, Quantity Desired, and Purchase Order Number.
THANK YOU. YOUR BUSINESS IS APPRECIATED.