

## ANTIBODY SPECIFICATION SHEET

<b>PRODUCT:</b>	Monoclonal Benzoyl Ecgonine
<b>CATALOG NUMBER:</b>	P01-99-11M-P
<b>LOT NUMBER:</b>	S02
<b>FORMAT:</b>	Protein A Purified Material
<b>HOST:</b>	Murine
<b>ANTIGEN:</b>	Benzoyl Ecgonine-BSA
<b>ANTIBODY SPECIFICITY:</b>	Benzoyl Ecgonine and Cocaine
<b>CONCENTRATION:</b>	5.37 mg/ml by A280nm
<b>ISOTYPE:</b>	IgG-1
<b>BUFFER:</b>	0.05M Phosphate, 0.85% Sodium Chloride, pH 7.15
<b>PRESERVATIVE:</b>	0.05% Sodium Azide
<b>STORAGE:</b>	Short Term -- Refrigerated Long Term – Frozen

***MANUFACTURED BY BIOSTRIDE, INC.  
For Research and Manufacturer Use Only***

### **ORDERING INFORMATION**

Orders may be placed by phone, fax, or mail.  
Please provide Part Number, Description, Quantity Desired, and Purchase Order Number.  
THANK YOU. YOUR BUSINESS IS APPRECIATED.



**EastCoast Bio, Inc.**  
211 Wells Street  
P.O. Box 489  
North Berwick, ME 03906

P 207.676.7639  
F 207.676.7658  
E [info@eastcoastbio.com](mailto:info@eastcoastbio.com)  
[www.eastcoastbio.com](http://www.eastcoastbio.com)