

## ANTIBODY SPECIFICATION SHEET

<b>PRODUCT:</b>	Monoclonal THC
<b>PART NUMBER:</b>	P01-99-21M-P
<b>LOT NUMBER:</b>	S15
<b>FORMAT:</b>	Protein A Purified Material
<b>HOST:</b>	Murine
<b>ANTIGEN:</b>	Delta-8-THC-BSA
<b>ANTIBODY SPECIFICITY:</b>	THC and Metabolites
<b>CONCENTRATION:</b>	6.0 mg/mL by A280nm
<b>ISOTYPE:</b>	IgG-1
<b>BUFFER:</b>	PBS-Azide
<b>PRESERVATIVE:</b>	0.05% Sodium Azide
<b>STORAGE:</b>	Short Term -- Refrigerated Long Term -- Frozen
<b>CROSS REACTIVITY:</b>	No cross reactivity with: <ul style="list-style-type: none"><li>• Opiates</li><li>• Cocaine Metabolite</li><li>• Amphetamines</li><li>• Phencyclidine</li></ul>

### ***MANUFACTURED BY BIOSTRIDE, INC.***

*This product is intended for research and manufacturing uses only. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.*

### **ORDERING INFORMATION**

Orders may be placed by phone, fax, or mail.  
Please provide Part Number, Description, Quantity Desired, and Purchase Order Number.  
THANK YOU. YOUR BUSINESS IS APPRECIATED.