

**ANTIBODY SPECIFICATION SHEET**

**PRODUCT:** Monoclonal THC  
**CATALOG NUMBER:** P01-99-22M-P  
**LOT NUMBER:** S11  
**FORMAT:** Protein A Purified Material  
**HOST:** Murine  
**ANTIGEN:** THC-BSA  
**ANTIBODY SPECIFICITY:** THC and Metabolite  
**CONCENTRATION:** 7.50 mg/ml  
by A280nm  
**ISOTYPE:** IgG-1  
**BUFFER:** PBS, pH 7.2  
**PRESERVATIVE:** 0.05% Sodium Azide  
**STORAGE:** Short Term -- Refrigerated  
Long Term -- Frozen

***MANUFACTURED BY BIOSTRIDE, INC.***

*This product is intended for research and manufacturing uses only. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations*

**ORDERING INFORMATION**

Orders may be placed by phone, fax, or mail.

Please provide Part Number, Description, Quantity Desired, and Purchase Order Number.  
THANK YOU. YOUR BUSINESS IS APPRECIATED.