

ANTIBODY SPECIFICATION SHEET

PRODUCT:	Polyclonal Fentanyl
CATALOG NUMBER:	P01-99-51R-IF
LOT NUMBER:	S05
FORMAT:	Purified Ig Fraction
HOST:	Rabbit
ANTIGEN:	Fentanyl-BgG (linked off the aliphatic group, ethyl)
ANTIBODY SPECIFICITY:	Fentanyl
CONCENTRATION:	3.49 mg/ml by A280nm
ISOTYPE:	N/A
BUFFER:	0.05 M Phosphate 0.85% Sodium Chloride, pH 7.15
PRESERVATIVE:	0.05% Sodium Azide
STORAGE:	Short Term -- Refrigerated Long Term – Frozen
COMMENTS:	Fentanyl linked off the aliphatic group, ethyl

***MANUFACTURED BY BIOSTRIDE, INC.
For Research and Manufacturing Use Only***

ORDERING INFORMATION

Orders may be placed by phone, fax, or mail.
Please provide Part Number, Description, Quantity Desired, and Purchase Order Number.
THANK YOU. YOUR BUSINESS IS APPRECIATED.



EastCoast Bio, Inc.
211 Wells Street
P.O. Box 489
North Berwick, ME 03906

P 207.676.7639
F 207.676.7658
E info@eastcoastbio.com
www.eastcoastbio.com