

**ANTIBODY SPECIFICATION SHEET**

|                              |   |
|------------------------------|---|
| <b>PRODUCT:</b>              | Monoclonal Phencyclidine  |
| <b>CATALOG NUMBER:</b>       | P01-99-65M-P  |
| <b>LOT NUMBER:</b>           | P04   |
| <b>FORMAT:</b>               | Protein A Purified Material   |
| <b>HOST:</b>                 | Murine  |
| <b>ANTIGEN:</b>              | Phencyclidine-BgG   |
| <b>ANTIBODY SPECIFICITY:</b> | Phencyclidine   |
| <b>CONCENTRATION:</b>        | 2.44 mg/ml  |
| <b>ISOTYPE:</b>              | IgG   |
| <b>BUFFER:</b>               | 0.05M Phosphate<br>0.85% Sodium Chloride, pH 7.15   |
| <b>PRESERVATIVE:</b>         | 0.05% Sodium Azide  |
| <b>STORAGE:</b>              | Short Term -- Refrigerated<br>Long Term – Frozen  |
| <b>CROSS REACTIVITY:</b>     | No apparent cross reactivity to:<br>Benzoyl Ecgonine<br>Opiates<br>Cannabinoids<br>Amphetamines |

***MANUFACTURED BY BIOSTRIDE, INC.  
for research and manufacturer use only***

**ORDERING INFORMATION**

Orders may be placed by phone, fax, or mail.  
Please provide Part Number, Description, Quantity Desired, and Purchase Order Number.  
THANK YOU. YOUR BUSINESS IS APPRECIATED.



**EastCoast Bio, Inc.**  
211 Wells Street  
P.O. Box 489  
North Berwick, ME 03906

P 207.676.7639  
F 207.676.7658  
E [info@eastcoastbio.com](mailto:info@eastcoastbio.com)  
[www.eastcoastbio.com](http://www.eastcoastbio.com)