

ANTIBODY SPECIFICATION SHEET

PRODUCT:	Monoclonal LSD
CATALOG NUMBER:	P01-99-66M-P
LOT NUMBER:	O23
FORMAT:	Neat Ascites and Purified
HOST:	Murine
ANTIGEN:	LSD- BgG
ANTIBODY SPECIFICITY:	LSD
CONCENTRATION:	1.39 mg/ml by A280nm
ISOTYPE:	IgG
BUFFER:	0.05M Potassium Phosphate, 0.85% NaCl, pH 7.15
PRESERVATIVE:	0.05% Sodium Azide
STORAGE:	Short Term -- Refrigerated Long Term -- Frozen

MANUFACTURED BY BIOSTRIDE, INC.

ORDERING INFORMATION

Orders may be placed by phone, fax, or mail.
Please provide Part Number, Description, Quantity Desired, and Purchase Order Number.
THANK YOU. YOUR BUSINESS IS APPRECIATED.



EastCoast Bio, Inc.
211 Wells Street
P.O. Box 489
North Berwick, ME 03906

P 207.676.7639
F 207.676.7658
E info@eastcoastbio.com
www.eastcoastbio.com



EastCoast Bio, Inc.
211 Wells Street
P.O. Box 489
North Berwick, ME 03906

P 207.676.7639
F 207.676.7658
E info@eastcoastbio.com
www.eastcoastbio.com