

ANTIBODY SPECIFICATION SHEET

PRODUCT: Monoclonal LSD

CATALOG NUMBER: P01-99-68M-P

LOT NUMBER: G48

FORMAT: Protein A Purified Material

HOST: Murine

ANTIGEN: LSD- BgG

ANTIBODY SPECIFICITY: LSD

CONCENTRATION: 2.62 mg/ml
by A280nm

ISOTYPE: IgG

BUFFER: 0.05M Phosphate
0.85% NaCl, pH 7.15

PRESERVATIVE: 0.05% Sodium Azide

STORAGE: Short Term -- Refrigerated
Long Term -- Frozen

MANUFACTURED BY BIOSTRIDE, INC.

ORDERING INFORMATION

Orders may be placed by phone, fax, or mail.
Please provide Part Number, Description, Quantity Desired, and Purchase Order Number.
THANK YOU. YOUR BUSINESS IS APPRECIATED.