

ANTIBODY SPECIFICATION SHEET

PRODUCT:	Monoclonal LSD
CATALOG NUMBER:	P01-99-69M-P
LOT NUMBER:	L01
FORMAT:	Protein A Purified Material
HOST:	Murine
ANTIGEN:	LSD- BgG
ANTIBODY SPECIFICITY:	LSD
CONCENTRATION:	2.20 mg/ml by A280nm
ISOTYPE:	IgG
BUFFER:	0.05M Phosphate 0.85% Sodium Chloride, pH 7.15
PRESERVATIVE:	0.05% Sodium Azide
STORAGE:	Short Term -- Refrigerated Long Term -- Frozen

MANUFACTURED BY BIOSTRIDE, INC.

ORDERING INFORMATION

Orders may be placed by phone, fax, or mail.
Please provide Part Number, Description, Quantity Desired, and Purchase Order Number.
THANK YOU. YOUR BUSINESS IS APPRECIATED.



EastCoast Bio, Inc.
211 Wells Street
P.O. Box 489
North Berwick, ME 03906

P 207.676.7639
F 207.676.7658
E info@eastcoastbio.com
www.eastcoastbio.com