

ANTIBODY SPECIFICATION SHEET

PRODUCT: Polyclonal Benzodiazepine
CATALOG NUMBER: P01-99-92S-IF
LOT NUMBER: 9/02
FORMAT: DEAE Purified Ig Fraction
HOST: Sheep
ANTIGEN: Oxazepam-BgG
ANTIBODY SPECIFICITY: Benzodiazepines
CONCENTRATION: 3.72 mg/mL
 By A280nm
ISOTYPE: N/A
BUFFER: 0.05 M. Phosphate
 0.85% Sodium Chloride, pH 7.15
PRESERVATIVE: 0.05% Sodium Azide
STORAGE: Short Term -- Refrigerated
 Long Term -- Frozen

CROSS REACTIVITY:	<u>ANALYTE</u>	<u>D50</u>	<u>AT 0.3ug/mL</u>
	Alprazolam	>100%	>100%
	Bromazepam	>100%	>100%
	Chlorazepate	>100%	>100%
	Chlordiazepoxide	12%	14%
	Clobazam	>100%	60%
	Clonazepam	50%	20%
	Delorazepam	>100%	69%
	D-methyldiazepam	>100%	>100%
	Diazepam	>100%	>100%
	Estazolam	>100%	>100%
	Flunitrazepam	>100%	50%
	Flurazepam	>100%	>100%
	Lorazepam	50%	21%
	Lormetazepan	>100%	40%
	Medazepam	>100%	>100%
	Nitrazepam	>100%	>100%
	Oxazolam	2%	1%
	Pramazepam	>100%	>100%
	Temazepam	>100%	>100%
	Triazolam	>100%	>100%

MANUFACTURED BY BIOSTRIDE, INC.

ORDERING INFORMATION

Orders may be placed by phone, fax, or mail.

Please provide Part Number, Description, Quantity Desired, and Purchase Order Number.

THANK YOU. YOUR BUSINESS IS APPRECIATED.