

ANTIBODY SPECIFICATION SHEET

PRODUCT: Polyclonal Benzodiazepine
CATALOG NUMBER: P01-99-93S-NS
LOT NUMBER: I32
FORMAT: Neat Serum
HOST: Sheep
ANTIGEN: Oxazepam-BgG
ANTIBODY SPECIFICITY: Benzodiazepines
ISOTYPE: N/A
STORAGE: Short Term -- Refrigerated
 Long Term -- Frozen

CROSS REACTIVITY:	<u>ANALYTE</u>	<u>D50</u>	<u>AT 0.3ug/mL</u>
	Alprazolam	>100%	>100%
	Bromazepam	>100%	>100%
	Chlorazepate	>100%	>100%
	Chlordiazepoxide	>100%	13%
	Clobazam	>100%	>100%
	Clonazepam	>100%	7%
	Delorazepam	>100%	15%
	D-methyl diazepam	>100%	>100%
	Diazepam	>100%	>100%
	Estazolam	>100%	>100%
	Flunitrazepam	>100%	50%
	Flurazepam	>100%	>100%
	Lorazepam	21%	4%
	Lormetazepam	>100%	7%
	Medazepam	>100%	>100%
	Nitrazepam	>100%	>100%
	Oxazolam	9%	4%
	Pramazepam	>100%	>100%
	Temazepam	>100%	>100%
	Triazolam	>100%	43%

MANUFACTURED BY BIOSTRIDE, INC.

ORDERING INFORMATION

Orders may be placed by phone, fax, or mail.

Please provide Part Number, Description, Quantity Desired, and Purchase Order Number.

THANK YOU. YOUR BUSINESS IS APPRECIATED.