

**ANTIBODY SPECIFICATION SHEET**

**PRODUCT:** Monoclonal Methadone

**CATALOG NUMBER:** P11-99-08M-P

**LOT NUMBER:** S33

**FORMAT:** Protein A Purified Material

**HOST:** Murine

**ANTIGEN:** Methadone-BSA

**ANTIBODY SPECIFICITY:** Methadone

**CONCENTRATION:** 7.43 mg/mL  
by A280nm

**ISOTYPE:** IgG

**BUFFER:** 0.05M Potassium Phosphate  
0.85% NaCl, pH 7.15

**PRESERVATIVE:** 0.05% Sodium Azide

**STORAGE:** Short Term -- Refrigerated  
Long Term -- Frozen

**CROSS REACTIVITY:** D50 in ELISA : Methadol ~ 2.9%

***MANUFACTURED BY BIOSTRIDE, INC.***

*This product is intended for research and manufacturing uses only. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.*

**ORDERING INFORMATION**

Orders may be placed by phone, fax, or mail.  
Please provide Part Number, Description, Quantity Desired, and Purchase Order Number.  
THANK YOU. YOUR BUSINESS IS APPRECIATED.