

**POLYCLONAL ANTIBODY DATA SHEET**

<b>PRODUCT:</b>	Anti-AGE (Advanced Glycosylation Endproduct)
<b>CATALOG NUMBER:</b>	PP102IF
<b>LOT NUMBER:</b>	R1891
<b>FORMAT:</b>	Ig Fraction, Liquid
<b>PURIFICATION:</b>	Protein A purified, 0.2 µm filtered
<b>PURITY:</b>	≥95% based on SDS-PAGE
<b>HOST:</b>	Rabbit
<b>CONCENTRATION:</b>	5.01 mg/ml (A280)
<b>IMMUNOGEN:</b>	Bovine Serum Albumin/Human Serum Albumin – Advanced Glycosylation Endproduct
<b>BUFFER:</b>	5 mM phosphate, 0.85% NaCl, @ pH 7.15
<b>CROSS REACTIVITY:</b>	BSA-AGE; HSA-AGE(< than 1%)
<b>PRESERVATIVES:</b>	0.05% Sodium Azide
<b>STORAGE:</b>	Short term: 4°C Long term: -20°C <i>Avoid repeated freeze-thaw cycles</i>

*This product is intended for research and manufacturing uses only. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.*

**ORDERING INFORMATION**

Orders may be placed by phone, fax, or mail.  
Please provide Part Number, Description, Quantity Desired, and Purchase Order Number.  
**THANK YOU. YOUR BUSINESS IS APPRECIATED.**



**EastCoast Bio, Inc.**  
211 Wells Street  
P.O. Box 489  
North Berwick, ME 03906

P 207.676.7639  
F 207.676.7658  
E [info@eastcoastbio.com](mailto:info@eastcoastbio.com)  
W [www.eastcoastbio.com](http://www.eastcoastbio.com)