

POLYCLONAL ANTIBODY DATA SHEET

PRODUCT: Anti-Human TSH

CATALOG NUMBER: PP135NS

LOT NUMBER: H2872

HOST: Goat

IMMUNOGEN: TSH

FORMAT: Neat Serum, liquid

PRESERVATIVES: None

STORAGE: Short term 4°C
Long term -20°C

This product is being sold for research purposes only. It is not to be used in humans or animals. There is no express or implied warranty. No liability is assumed for fitness of purpose and merchantability, or direct and consequential damages. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations. 10/03